

Anti-MSLN antibody (296-606) (STJ116034)

STJ116034

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IHC-P:1:50-1:200 IF/ICC: 1:100-1:500 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	10232
Gene Symbol	MSLN
Uniprot ID	MSLN_HUMAN
Immunogen	
Immunogen Region	296-606
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 296-606 of human Mesothelin (NP_037536.2).
Immunogen Sequence	EVEKTACPSGKKAREIDESL IFYKKWELEACVDAALLATQ MDRVNAIPFTYEQLDVLKHK LDELYPQGYPESVIQHLGYL FLKMSPEDIRKWNVTSLLETL KALLEVNKGHEMSPQAPRRP LPQVATLIDRFVKGRGQLDK DTLDTLTAFYPGYLCSLSPE ELSSVPPSSIWAVRPQDLDT CDPRQLDVLYPKARLAFQNM NGSEYFVKIQSFLGGAPTED LKALSQQNVSMDLATFMKL




Western blot analysis of extracts of OVCAR3 cells, using Mesothelin antibody (STJ116034) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS00855) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry analysis of paraffin-embedded human lung cancer (200X) using Mesothelin rabbit polyclonal antibody (STJ116034) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human lung cancer (400X) using Mesothelin rabbit polyclonal antibody (STJ116034) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded Human tonsil (200X) using Mesothelin rabbit polyclonal antibody (STJ116034) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081