

Anti-CACNG2 antibody (204-323) (STJ116036)

STJ116036

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10369
Gene Symbol	CACNG2
Uniprot ID	CCG2_HUMAN
Immunogen	
Immunogen Region	204-323
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 204-323 of human CACNG2 (NP_006069.1).
Immunogen Sequence	RHKQLRATARATDYLAQSAI TRIPSYRYRYQRRSRSSRS TEPHSRDASPVGIKGFNTL PSTEISMYTLSDPLKAATT PTATYNSDRDNSFLQVHNCI QKENKDSLHSNTANRRTPV