

Anti-ERCC6L antibody (300-400) (STJ116072)

STJ116072

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54821
Gene Symbol	ERCC6L
Uniprot ID	ERCC6L_HUMAN
Immunogen	
Immunogen Region	300-400
Specificity	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ERCC6L (NP_060139.2).
Immunogen Sequence	KDATPGEKALGFKISENLMA IIKPYFLRRTKEDVQKKSS NPEARLNEKNPDVDAICEMP SLSRKNDLIWIRLVPLQEE IYRKFSLDHIKELLMETRS P