

Anti-SETD5 antibody (520-740) (STJ116075)

STJ116075

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55209
Gene Symbol	SETD5
Uniprot ID	SETD5_HUMAN
Immunogen	
Immunogen Region	520-740
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 520-740 of human SETD5 (NP_001073986.1).
Immunogen Sequence	HAFENLEKRRKKRRDQPLEQS NSDVEITTTTSETPVGEETK TEAPESEVSNVSNVTIPST PQSVGVNTRRSSQAGDIAAE KLVKPPPAKPSRPRPKSRI SRYRTSSAQLRKQKQANAQ QAELSQAALKEEGGNSLVTP TEAGSLDSSGENRPLTGSDP TVVSITGSHVNRAASKYPKT KKYLVTWELNDKAEKQECVP ECPRLITTDPTVLATTLNML P



Immunofluorescence analysis of MCF-7 cells using SETD5 antibody (STJ116075). Blue: DAPI for nuclear staining.