

Anti-GSTO2 antibody (124-243) (STJ116091)

STJ116091

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/CC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	119391
Gene Symbol	GSTO2
Uniprot ID	GSTO2_HUMAN
Immunogen	
Immunogen Region	124-243
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 124-243 of human GSTO2 (NP_899062.1).
Immunogen Sequence	PHLTKECLVALRCGRECTNL KAALRQEFNSNLEEILEYQNT TFFGGTCISMIDYLLWPWFE RLDVYGILDCVSHTPALRLW ISAMKWDPTVCALLMDKSIF QGFLNLYFQNNPNAFDGLC



Immunofluorescence analysis of U2OS cells using GSTO2 antibody (STJ116091). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081