

Anti-SERPINA12 antibody (20-320) (STJ116095)

STJ116095

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	145264
Gene Symbol	SERPINA12
Uniprot ID	SPA12_HUMAN
Immunogen	
Immunogen Region	20-320
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human SERPINA12 (NP_776249.1).
Immunogen Sequence	LLKPSFSRNYKALSEVQGW KQRMAAKELARQNMDLGFKL LKKLAFYNPGRNIFLSPLSI STAFSMLCLGAQDSTLDEIK QGFNFRKMPEKDLHEGFHYI IHELTQKTQDLKLSIGNTLF IDQRLQPQRKFLEDAKNFYS AETILTQNFQNLQMAQKQIND FISQKTHGKINNLIENIDPG TVMLLANYIFFRARWKHEFD PNVTKEDFFLEKNSSVKVP MMFRSGIYQVGYDDKLSCT