

Anti-KBTBD3 antibody (410-607) (STJ116168)

STJ116168

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	69149
Gene Symbol	Kbtbd3
Uniprot ID	KBTB3_MOUSE
Immunogen	
Immunogen Region	410-607
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 410-607 of mouse KBTBD3 (NP_081238.2).
Immunogen Sequence	GSQDIKSLLDVESYNPLSRE WTSVSPLPRGIYYPEASACQ NIIYVLGSEVEIADAFNPSL DCFFKYNNATTDQWSELVAEF GQFFHATLIKAVPVNCTLYI CDLSTYKVVYFCPDTVCVWKG EGSFECAGFNAGAIGIEDKI YILGGDYAPDEITDEVQVYH SSRSEWEEVSPMPRALTEFY CQVIQFNKYRDPWYSNHF



Western blot analysis of extracts of PC-3 cells, using KBTBD3 antibody (STJ116168). Secondary antibody: HRP-Goat Anti-Rabbit IgG (H+L) (STJSD00856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081