

Anti-REST antibody (760-1060) (STJ116170)

STJ116170

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5978
Gene Symbol	REST
Uniprot ID	REST_HUMAN
Immunogen	
Immunogen Region	760-1060
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 760-1060 of human REST/NRSF (NP_001180437.1).
Immunogen Sequence	PVKIELSPPIEVQKEPVQM ELSPPMGVVQKEPAQREPPP PREPPLHMEPISKKPPLRKD KKEKSNMQSERARKEQVLIE VGLVPVKDSSLWKESVSTED LSPSPPLPKENLREEASGD QKLLNTGEGNKEAPLQKVGA EEADESLPGLAANINESTHI SSSGQNLNTPEGETLNGKHQ TDSIVCEMKMDTDQNTRENL TGINSTVEEPVSPMLPPSAV EEREAVSKTALASPPATMA



Immunofluorescence analysis of MCF-7 cells using REST/NRSF Rabbit polyclonal antibody (STJ116170). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081