

Anti-LEPREL2 antibody (487-736) (STJ116174)
STJ116174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10536
Gene Symbol	P3H3
Uniprot ID	P3H3_HUMAN
Immunogen Region	487-736
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 487-736 of human LEPREL2 (NP_055077.2).
Immunogen Sequence	DAAGAGARSGYRGRSPHTP HERFEGTLVLKAAQLARAGT VGSQGAKLLLEVSRVRLT QAYFSPERPLHLSFTHLVCRA SAIEGEQEQRMDLSHPVHAD NCVLDPDTGECWREPPAYTY RDYSGLLYLNDDFQGGDLFF TEPNALTVTARVRPRCGRLV AFSSGVENPHGWAVTRGR CALALWHTWAPEHREQEWIE AKELLQESQEEEEEEEMP SKDPSPEPPSRRHQRVQDK