

Anti-DACH1 antibody (450-570) (STJ116217)

STJ116217

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1602
Gene Symbol	DACH1
Uniprot ID	DACH1_HUMAN
Immunogen	
Immunogen Region	450-570
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 450-570 of human DACH1 (NP_542937.2).
Immunogen Sequence	EGRRPGSHPSHRSSSVSSS PARTESSDRIPVHQGLSM NQMLMGLSPNVLPGPKEGDL AGHDMGHESKRMHIEKDETP LSTPTARDSLDKLSLTGHGQ PLPPGFPPFLFPDGLSSIE T