

Anti-IWS1 antibody (570-819) (STJ116263)

STJ116263

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55677
Gene Symbol	IWS1
Uniprot ID	IWS1_HUMAN
Immunogen	
Immunogen Region	570-819
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 570-819 of human IWS1 (NP_060439.2).
Immunogen Sequence	AAEEDRQLNNQKKPALKKLT LLPAVVMHLKKQDLKETFD SGVMSAIKEWLSPLPDRSLP ALKIREELLKILQELPSVSQ ETLKHSGIGRAVMYLYKHPK ESRSNKDMAGKLINEWSRPI FGLTSNYKGMTREEREQRDL EQMPQRRRMNSTGGQTPrRD LEKVLTGEEKALRPGDPGFC ARARVPMPSNKDYVVRPKWN VEMESSRFQATSKKGISRLD KQMRKFTDIRKSRSAHAV