

Anti-EIF5B antibody (901-1220) (STJ116285)

STJ116285

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9669
Gene Symbol	EIF5B
Uniprot ID	IF2P_HUMAN
Immunogen Region	901-1220
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1220 of human EIF5B (NP_056988.3).
Immunogen Sequence	IRGLLLPPPMKELRVKNQYE KHKEVEAAQGVKILGKDLEK TLAGLPLLVAYPEDEIPVLK DELIHELKQTLNAIKLEEKG VYVQASTLGSLEALLEFLKT SEVPYAGINIGPVHKKDVMK ASVMLEHDPQYAVILAFDVR IERDAQEMADSLGVRIFSAE IIYHLDAFTKYRQDYKKQK QEEFKHIAVFPCKIKILPQY IFNSRDPVIMMGVTVEAGQVK QGTPMCVPSKNFVDIGIVT