

Anti-TRIM13 antibody (40-310) (STJ116298)
STJ116298

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10206
Gene Symbol	TRIM13
Uniprot ID	TRI13_HUMAN
Immunogen Region	40-310
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 40-310 of human TRIM13 (NP_001007279.1).
Immunogen Sequence	ILEGSRNSLWRPAPFKCPT CRKETSATGINSLQVNYSLK GIVEKYNKIKISPKMPVCKG HLGQPLNIFCLTDMQLICGI CATRGEHTKHVFCSDAY A QERDAFESLFQSFETWRRGD ALSRLDTLET SKRKSLLT KDSDKVKEFFEKLQHTLDQK KNEILSDFETMKLAVMQAYD PEINKLNTILQEQRMAFNIA EAFKDVSEPIVFLQQMQEFR EKIKVIKETPLPPSNLPAS