

Anti-SPAG5 antibody (1-300) (STJ116304)

STJ116304

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:20-1:50
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10615
Gene Symbol	SPAG5
Uniprot ID	SPAG5_HUMAN
Immunogen	
Immunogen Region	1-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SPAG5 (NP_006452.3).
Immunogen Sequence	MWRVKKLSLSLSPSPQTGKP SMRTPLRELTLPGLTNSG KRSPACSSLTPSLCKLGLQE GSNNSSPVDFVNNKRTDLSS EHFSHSSKWLETQHESEDEQ PLDPIPIQSSTPKTSEEAVD PLGNMVKITIVLPSPLGQQ QDMIFEARLDTMAETNSISL NGPLRTDDLVRREEVAPCMGD RFSEVAAVSEKPIFQESPSH LLEESPPNPCSEQLHCSEKES LSSRTEAVREDLVPSESNA



Immunofluorescence analysis of HeLa cells using SPAG5 antibody (STJ116304). Blue: DAPI for nuclear staining.