

Anti-ITIH3 antibody (35-180) (STJ116359)

STJ116359

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3699
Gene Symbol	ITIH3
Uniprot ID	ITIH3_HUMAN
Immunogen	
Immunogen Region	35-180
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 35-180 of human ITIH3 (NP_002208.3).
Immunogen Sequence	LPEGVANGIEVYSTKINSKV TSRFAHNVVTMRAVNRADTA KEVSFDVELPKTAFITNFTL TIDGVTPGNVKEEVAKKQ YEKAVSQGKTAGLVKASGRK LEKFTVSVNVAAGSKVTFEL TYEELLKRHKGYEMYLVKVQ PKQLVK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081