

Anti-ARSB antibody (404-533) (STJ116444)

STJ116444

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	411
Gene Symbol	ARSB
Uniprot ID	ARSB_HUMAN
Immunogen	
Immunogen Region	404-533
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 404-533 of human ARSB (NP_000037.2).
Immunogen Sequence	PCPRNSMAPAKDDSSLPEYS AFNTSVHAAIRHGNWKLITG YPGCGYWFPPPSQYNVSEIP SSDPPTKTLWLFDIDRDPEE RHDLSREYPHIVTKLLSRLQ FYHKHSVPVYFPAQDPRCDP KATGWGWPWM



Western blot analysis of lysates from U-87MG cells, using ARSB Rabbit polyclonal antibody (STJ116444) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081