

Anti-CAND1 antibody (750-980) (STJ116499)

STJ116499

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	55832
Gene Symbol	CAND1
Uniprot ID	CAND1_HUMAN
Immunogen	
Immunogen Region	750-980
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 750-980 of human CAND1 (NP_060918.2).
Immunogen Sequence	GALSAMLDFFQALVVTGNN LGYMDLLRMLTGPVYSQSTA LTHKQSYYSIAKCVAAALTRA CPKEGPAVVGQFIQDVKNRSR STDSIRLLALLSLGEVGHII DLSGQLELKSVILEAFSSPS EEVKSAAAYALGSISVGNLP EYLPFVLQETISQPKRQYLL LHSLKEISSASVGLKPYV ENIWALLLKHCCEAEEGTRN VVAECLGKLTLDPETLLPR LKGYLISGSSY




Immunohistochemistry analysis of CAND1 in paraffin-embedded human uterine cancer using CAND1 Rabbit polyclonal antibody (STJ116499) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of CAND1 in paraffin-embedded human colon using CAND1 Rabbit polyclonal antibody (STJ116499) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of CAND1 in paraffin-embedded human lung cancer using CAND1 Rabbit polyclonal antibody (STJ116499) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of CAND1 in paraffin-embedded human lung using CAND1 Rabbit polyclonal antibody (STJ116499) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081