

Anti-RHBDD1 antibody (206-315) (STJ116576)

STJ116576

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84236
Gene Symbol	RHBDD1
Uniprot ID	RHBL4_HUMAN
Immunogen	
Immunogen Region	206-315
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 206-315 of human RHBDD1 (NP_115652.2).
Immunogen Sequence	TQGPKKIMEACAGGFSSV GYPGRQYYFNSSGSGYQDY YPHGRPDHYEEAPRNYDTYT AGLSEEEQLERALQASLWDR GNTRNSPPPYGFHLSPEEMR RQRLHRFDSQ