

Anti-ZNF286A antibody (120-240) (STJ116625)

STJ116625

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57335
Gene Symbol	ZNF286A
Uniprot ID	Z286A_HUMAN
Immunogen	
Immunogen Region	120-240
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 120-240 of human ZNF286A (NP_065703.1).
Immunogen Sequence	KDSTSVQDFSKAESCKVAII DRLTRNSVYDSNLEAALECE NWLENQQGNQERHLREMFTH MNSLSEETDHHKHDVYWKSFN QKSVLITEDRVPKGSYAFHT LEKSLKQKSNLMKKQRTYKE K