

## Anti-SUMF1 antibody (34-374) (STJ116673)

STJ116673

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution</b>	WB:1:500-1:2000
<b>Range</b>	ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">285362</a>
<b>Gene Symbol</b>	<a href="#">SUMF1</a>
<b>Uniprot ID</b>	<a href="#">SUMF1_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	34-374
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 34-374 of human SUMF1 (NP_877437.2).
<b>Immunogen Sequence</b>	SQEAGTGAGAGSLAGSCGCG TPQRPGAHGSSAAAHYSRE ANAPGPVPGERQLAHSKMVP IPAGVFTMGTDDEPQIKODGE APARRVTIDAFYMDAYEVSN TEFEKFNSTGYLAEKFG DSFVFEGMLSEQVKTNIQQA VAAAPWWLPVKGANWRHPEG PDSTILHRPDHPVLHVSND AVAYCTWAGKRLPTEAEWEY SCRGGHNRLLFPWGNKLQPK GQHYANIWQGEFPVTNTGE