

Anti-PDIA5 antibody (22-140) (STJ116686)

STJ116686

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10954
Gene Symbol	PDIA5
Uniprot ID	PDIA5_HUMAN
Immunogen Region	22-140
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-140 of human PDIA5 (NP_006801.1).
Immunogen Sequence	SSAKVSSLIERISDPKDLKK LLRTRNNVLVLYSKSEVAE NHLRLLSTVAQAVKGGTIC WDCGDAESRKLCKMKVDL SPKDKKVELFHYQDGAFHTE YNRAVTFKSIVAFLKDPKG