

Anti-ASF1A antibody (100-204) (STJ116719)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25842
Gene Symbol	ASF1A
Uniprot ID	ASF1A_HUMAN
Immunogen	
Immunogen Region	100-204
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 100-204 of human ASF1A (NP_054753.1).
Immunogen Sequence	TYRGQEFIRVGYYVNNYEY TELRENPPVKPDFSKLQRNI LASNPRVTRFHINWEDNTEK LEDAESSNPNLQSLSTDAL PSASKGWSTSENSLNVMLCS HMDCM