

Anti-DHRS11 antibody (1-260) (STJ116807)

STJ116807

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	79154
Gene Symbol	DHRS11
Uniprot ID	DHR11_HUMAN
Immunogen	
Immunogen Region	1-260
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DHRS11 (NP_077284.2).
Immunogen Sequence	MARPGMERWRDRLALVTGAS GGIGAARALVQQGLKVVG CARTVGNIEELAAECKSAGY PGTLPYRCDSLNEEDILSM FSAIRSQHSGVDICINNAGL ARPDTLLSGSTSGWKDMFNV NVLALSICTREAYQSMKERN VDDGHIININMSGHRVLPL SVTHFYSATKYAVTALTEGL RQELREAQTHIRATCISPGV VETQFAFKLHDKDPEKAAAT YEQMKCLKPEDVAEAVIYV