

Anti-CCT6B antibody (120-340) (STJ116824)
STJ116824

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:2000 IF/CC: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10693
Gene Symbol	CCT6B
Uniprot ID	TCPW_HUMAN
Immunogen	
Immunogen Region	120-340
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 120-340 of human CCT6B (NP_006575.2).
Immunogen Sequence	AEGFEAAKIKALEVLEEVKV TKEMKRKILLDVARTSLQTK VHAELADVLTEVVVDSVLAV RRPGYPIDLFMVEIMEMKHK LGTDTKLIQGLVLDHGARHP DMKKRVEDAFILCNVSLEY EKTEVNSGFFYKTAEEKEKL VKAERKFIEDRVQKIIDLKD KVCAQSNKGFFVINQKGIDP FSLDSLAKHGIVALRAKRR NMERLSLACGGMVAVNSFEDL T



Immunofluorescence analysis of L929 cells using CCT6B Rabbit polyclonal antibody (STJ116824) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of HeLa cells using CCT6B Rabbit polyclonal antibody (STJ116824) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using CCT6B Rabbit polyclonal antibody (STJ116824) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081