

Anti-IFT57 antibody (170-429) (STJ116858)

STJ116858

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55081
Gene Symbol	IFT57
Uniprot ID	IFT57_HUMAN
Immunogen	
Immunogen Region	170-429
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 170-429 of human IFT57 (NP_060480.1).
Immunogen Sequence	VEELEESVAEDDAELTLNK VDEEFVEEETDNEENFIDLN VLKAQTYHLDMMNETAKQEDI LESTTDAAEWSLEVERVLPQ LKVTIRTDNKDWRIHVDQMH QHRSGIESALKETKGFLDKL HNEITRTLKISSREKYINN QLENLVQEYRAAQAQLSEAK ERYQQGNGGVTERTRLLSEV MEELEKVKQEMEEKGSSMTD GAPLVKIKQSLTKLKQETVE MDIRIGIVEHTLLQSKLKE



Immunohistochemistry of paraffin-embedded human uterine cancer using IFT57 antibody (STJ116858) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat ovary using IFT57 antibody (STJ116858) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081