

Anti-PAEP antibody (19-180) (STJ116957)

STJ116957

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5047
Gene Symbol	PAEP
Uniprot ID	PAEP_HUMAN
Immunogen	
Immunogen Region	19-180
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-180 of human PAEP (NP_002562.2).
Immunogen Sequence	MDIPQTKQDLELPKLAGTWH SMAMATNNISLMATLKAPLR VHITSLPTPEDNLEIVLHR WENNSCVEKKVLGEKTENPK KFKINYTVANEATLLD TDYD NFLFLCLQD TTTPIQSMM CQ YLARVLVEDDEIMQGFIRAF RPLPRHLWYLLDLKQMEEPC RF