

Anti-TFDP1 antibody (1-290) (STJ116989)

STJ116989

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7027
Gene Symbol	TFDP1
Uniprot ID	TFDP1_HUMAN
Immunogen	
Immunogen Region	1-290
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human DP1/DP1/TFDP1 (NP_009042.1).
Immunogen Sequence	MAKDAGLIEANGELKVFIDQ NLSPGKGVVSLVAVHPSTVN PLGKQLLPKTFGQSNVNIAQ QVWIGTPQRPAAASNTLVVGS PHTPSTHFASQNPQSDSSPW SAGKRNKGEKNGKGLRHFS MKVCEKVQRKGTTTSYNEVAD ELVAEFSAADNHILPNESAY DQKNIRRRVYDALNVLAMN IISKEKKEIKWIGLPTNSAQ ECQNLEVERQRRLERIKQKQ SQLQELIQQIAFKNLVQR



Immunohistochemistry analysis of paraffin-embedded rat testis using DP1/DP1/TFDP1 antibody (STJ116989) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.