

Anti-KCNAB3 antibody (285-404) (STJ117021)

STJ117021

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9196
Gene Symbol	KCNAB3
Uniprot ID	KCAB3_HUMAN
Immunogen	
Immunogen Region	285-404
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 285-404 of human KCNAB3 (NP_004723.2).
Immunogen Sequence	YPLACGLITSKYDGRVPDTC RASIKGYQWLKDKVQSEDGK KQAKVMDLLPVAHQQLGCTV AQLAIWCLRSEGVSSVLLG VSSAEQLIEHLGALQVLSQL TPQTVMEIDGLLGKPKPHSKK



Western blot analysis of various lysates using KCNAB3 Rabbit polyclonal antibody (STJ117021) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081