

Anti-DCTN2 antibody (280-401) (STJ117041)

STJ117041

GENERAL INFORMATION


| | |
|--------------------------|--------------------|
| Product Type | Primary antibodies |
| Short Description | |
| Applications | WB/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | Lot specific |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 10540 |
| Gene Symbol | DCTN2 |
| Uniprot ID | DCTN2_HUMAN |
| Immunogen | |
| Immunogen Region | 280-401 |
| Specificity | Recombinant fusion protein containing a sequence corresponding to amino acids 280-401 of human DCTN2 (NP_006391.1). VLDQVEARLQSVLGKVNEIA KHKASVEDADTQSKVHQLYE TIQRWSPIASTLPVLVQRLV TIKQLHEQAMQFGQLLTHLD |
| Immunogen Sequence | TTQQMIA NSLKDNTLLTQV QTTMRENLATVEGNFASIDE RM |



Western blot analysis of various lysates, using DCTN2 Rabbit polyclonal antibody (STJ117041) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081