

Anti-ITGB3BP antibody (1-111) (STJ117059)
STJ117059

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	23421
Gene Symbol	ITGB3BP
Uniprot ID	CENPR_HUMAN
Immunogen	
Immunogen Region	1-111
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human ITGB3BP (NP_055103.3).
Immunogen Sequence	MPVKRSLKLDGLLENSFDP SKITRKKSIVITYSPTTGTCQ MSLFASPTSSEEQKHRNGLS NEKRKKLNHPSLTESKESTT KDNEFMMLLSKVEKLSEEI MEIMQNLSIQ




Immunohistochemistry analysis of ITGB3BP in paraffin-embedded mouse brain using ITGB3BP Rabbit polyclonal antibody (STJ117059) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of ITGB3BP in paraffin-embedded human gastric cancer using ITGB3BP Rabbit polyclonal antibody (STJ117059) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of ITGB3BP in paraffin-embedded human appendix using ITGB3BP Rabbit polyclonal antibody (STJ117059) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of ITGB3BP in paraffin-embedded human liver using ITGB3BP Rabbit polyclonal antibody (STJ117059) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081