

Anti-FBXO42 antibody (468-717) (STJ117098)

STJ117098

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	54455
Gene Symbol	FBXO42
Uniprot ID	FBX42_HUMAN
Immunogen	
Immunogen Region	468-717
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 468-717 of human FBXO42 (NP_061867.1).
Immunogen Sequence	PSAPEGYDLKIGLSLAPRRG SLPDQKDLRLGSIDLNWDLK PASSSNPMDGMDNRTVGGSM RHPPEQTNGVHTPPHVASAL AGAVSPGALRRSLEAIKAMS SKGPSAASALSPPLGSSPGS PGSQSLSSGETVPIPRPGPA QGDGHSPLPIARRLGHHPQ SLNVGKPLYQSMNCKPMQMY VLDIKDTEKGRVKWKVFNS SSVVGPPETSLHTVVQGRGE LIIFGGLMDKKQNVKYYPK



Western blot analysis of extracts of mouse lung, using FBXO42 antibody (STJ117098) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081