

## Anti-ZNF415 antibody (100-260) (STJ117113)

STJ117113

### GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IHC-P/ELISA     |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Clone ID</b>            |                                                                                                                                        |
| <b>Concentration</b>       | Lot specific                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Purification</b>        | Affinity purification                                                                                                                  |
| <b>Dilution</b>            | WB:1:500-1:2000                                                                                                                        |
| <b>Range</b>               | IHC-P:1:50-1:200                                                                                                                       |
|                            | ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.                                                                                     |
| <b>Isotype</b>             | IgG                                                                                                                                    |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                         |

### TARGET INFORMATION


|                           |                                                                                                                                                                        |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 55786                                                                                                                                                                  |
| <b>Gene Symbol</b>        | ZNF415                                                                                                                                                                 |
| <b>Uniprot ID</b>         | ZN415_HUMAN                                                                                                                                                            |
| <b>Immunogen</b>          |                                                                                                                                                                        |
| <b>Immunogen Region</b>   | 100-260                                                                                                                                                                |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human ZNF415 (NP_060825.2).                                                   |
| <b>Immunogen Sequence</b> | WRDDEKNCNKVTTAPKENLT CRRDQRDRRGIGNKSIKHL GLSFLPHPELQQFQAEGKI YECNHVEKSVNHGSSVSPQ IISSTIKTHVSNKYGTDFIC SLLLTQEQKSCIREKPYRYI ECDKALNHGSHMTVRQVSHS GEKGYKCDLCKGVFSQKSNL A |



Immunohistochemistry of paraffin-embedded rat brain using ZNF415 antibody (STJ117113) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using ZNF415 antibody (STJ117113) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using ZNF415 antibody (STJ117113) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081