

Anti-AMICA1 antibody (285-384) (STJ117157)

STJ117157

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	120425
Gene Symbol	JAML
Uniprot ID	JAML_HUMAN
Immunogen	
Immunogen Region	285-384
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 285-384 of human AMICA1 (NP_694938.2).
Immunogen Sequence	IVKKTCGNKSSVNSTVLVKN TKKTNPETKEKPCFERCEG EKHIYSPIIVREVIEEEPS EKSEATYMTMHPVWPSLRSD RNNSLEKKSGGMPKTQQAF