

Anti-ADAM8 antibody (17-300) (STJ117216)

STJ117216

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	101
Gene Symbol	ADAM8
Uniprot ID	ADAM8_HUMAN
Immunogen	
Immunogen Region	17-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human ADAM8 (NP_001100.3).
Immunogen Sequence	IAPSRPWALMEQYEVLPWR LPGPRVRRALPSHLGLHPER VSYVLGATGHNFTLHLRKNR DLLGSGYTETYTAANGSEVT EQPRGQDHCIFYQGHVEGYPD SAASLSTCAGLRGFFQVGSD LHLIEPLDEGGEGGRHAVYQ AEHLLQTAGTCGVSDDSLGS LLGPRTAAVFRPRPGDSLPS RETRYVELYVVVDNAEFQML GSEAAVRHRVLEVNHVDKL YQKLNFRVVLVGLIEWNSQ