

Anti-ARSF antibody (1-200) (STJ117225)

STJ117225

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	416
Gene Symbol	ARSF
Uniprot ID	ARSF_HUMAN
Immunogen	
Immunogen Region	1-200
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ARSF (NP_004033.2).
Immunogen Sequence	MRPRRLVFMVSLVCALLNTC QAHRVHDDKPNIVLIMVDDL GIGDLGCYGNMTMRTPHIDR LAREGVRLTQHISAASLCSP SRSAFLTGRYPISGMVSSG NRRVIQNLAVPAGLPLNETT LAALLKKQGYSTGLIGKWHQ GLNCDSRSDQCHHPYNYGFD YYYGMPFTLVDSWPDPSRN TELAFESQLWLCVQLVAIAI



Immunofluorescence analysis of U2OS cells using ARSF antibody (STJ117225) at dilution of 1:100. Blue: DAPI for nuclear staining.