

Anti-CYP2A6 antibody (91-300) (STJ117238)

STJ117238

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1548
Gene Symbol	CYP2A6
Uniprot ID	CP2A6_HUMAN
Immunogen	
Immunogen Region	91-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 91-300 of human CYP2A6 (NP_000753.3).
Immunogen Sequence	LVDQAEFSGRGEQATFDWV FKGYGVSNGERAKQLRRF SIATLRDFGVGKRGIEERIQ EEAGFLIDALRGTTGGANIDP TFFLSRTVSNVISSIVFGDR FDYKDKEFLSLLRMMLGIFQ FTSTSTGQLYEMFSSVMKHL PGPQQQAFQLQGLEDFAK KVEHNQRTLDPNSPRDFDS FLIRMQEEKNPNTEFYLNK LVMTTNLNFI