

Anti-PSG2 antibody (200-300) (STJ117280)

STJ117280

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5670
Gene Symbol	PSG2
Uniprot ID	PSG2_HUMAN
Immunogen	
Immunogen Region	200-300
Specificity	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PSG2 (NP_112536.2).
Immunogen Sequence	RTLFLFGVTKYTAGPYECEI RNSGSASRSDPVTLLHGP DLPRIHPSYTNYRSGDNLYL SCFANSNPPAQYSWTINGKF QQSGQNLFIPIQITTKHSGLY V