

Anti-FZD4 antibody (37-222) (STJ117307)
STJ117307

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8322
Gene Symbol	FZD4
Uniprot ID	FZD4_HUMAN
Immunogen Region	37-222
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 37-222 of human FZD4 (NP_036325.2).
Immunogen Sequence	FGDEEERRCDPIRISM CQNL GYNVTKMPNLVGHELQTD AE LQLTFTPLIQYGCSSQLQF FLCSVYVPMCTEKINIPGP CGGMCLSVKRRCEPVLKEFG FAWPESLNC SKFPPQNDHNH MCMEGPGDEEVPLPHKTPIQ PGEECHSVGTNSDQYIWVKR SLNCVLKCGYDAGLYRS AK EFTDIW