

Anti-UBA3 antibody (234-463) (STJ117313)

STJ117313

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9039
Gene Symbol	UBA3
Uniprot ID	UBA3_HUMAN
Immunogen	
Immunogen Region	234-463
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 234-463 of human UBA3 (NP_003959.3).
Immunogen Sequence	FPMCTIASMPRLPEHCIEYV RMLQWPKEQPFGEVPLDGD DPEHIQWIFQKSLERASQYN IRGVTYRLTQGVVKRIIPAV ASTNAVIAAVCATEVFKIAT SAYIPLNNYLVFNDVDGLYT YTFEAERKENC PACSQLPON IQFSPSAKLQEVLDYLTNSA SLQMKSPAITATLEGKNRTL YLQSVTSIEERTRPNLSKTL KELGLVDGQELAVADVTPQ TVLFKLHFTS