

Anti-ZBTB20 antibody (100-400) (STJ117343)

STJ117343

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26137
Gene Symbol	ZBTB20
Uniprot ID	ZBT20_HUMAN
Immunogen	
Immunogen Region	100-400
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human ZBTB20 (NP_056457.3).
Immunogen Sequence	ILTAASILQIKTVIDECTRI VSQNVGDVFPGIQDSGQDTP RGTPESGTSGQSSDTESGYL QSHPOHSVDRIYSALYACSM QNGSGERSFYSGAVSHHET ALGLPRDHHMEDPSWITRIH ERSQQMERYLSTTPETTHCR KQPRPVRIQTLVGNIIHKQE MEDDYDYYGQQRVQILERNE SEECTEDTDQAEGTESEPKG ESFDSGVSSSIGTEPDSVEQ QFGPGAARDSQAEPQPEQ