

Anti-ZNF408 antibody (1-240) (STJ117384)

STJ117384

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79797
Gene Symbol	ZNF408
Uniprot ID	ZN408_HUMAN
Immunogen	
Immunogen Region	1-240
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ZNF408 (NP_079017.1).
Immunogen Sequence	MEEAEELLLEGKKALQLARE PRLGLDLGWNPSGEGCTQGL KDVPEPTRDILALKSLPRG LALGPSLAKEQRLGWVCVGD PLQPGLLWGPLEEESASKEK GEGVKPRQEENLSLGPWGDV CACEQSSGWTSLVQRGRLES EGNVAPVRISERLHLQVYQL VLPGSELLLWPQPSSEGPSL TQPGLDKEAAVAWTEVESA VQQEVASPGEDAAEPCIDPG SQSPSGIQAENMVSPGLKF



Immunofluorescence analysis of HeLa cells using ZNF408 antibody (STJ117384) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.