

# Anti-MYH2 antibody (621-635) (STJ117487)

STJ117487

## GENERAL INFORMATION




|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IF/ICC/ELISA    |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

## PRODUCT PROPERTIES


|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                          |
| <b>Clone ID</b>            |                                                                                                                                     |
| <b>Concentration</b>       | Lot specific                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                        |
| <b>Purification</b>        | Affinity purification                                                                                                               |
| <b>Dilution</b>            | WB:1:500-1:1000                                                                                                                     |
| <b>Range</b>               | IF/CC:1:50-1:200                                                                                                                    |
|                            | ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.                                                                                    |
| <b>Isotype</b>             | IgG                                                                                                                                 |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                      |

## TARGET INFORMATION

|                           |                                                                                                                                                                                                                                                                   |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 4620                                                                                                                                                                                                                                                              |
| <b>Gene Symbol</b>        | MYH2                                                                                                                                                                                                                                                              |
| <b>Uniprot ID</b>         | MYH2_HUMAN                                                                                                                                                                                                                                                        |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                                   |
| <b>Immunogen Region</b>   | 621-635                                                                                                                                                                                                                                                           |
| <b>Specificity</b>        | A synthetic peptide corresponding to a sequence within amino acids 600-700 of human MYH2 (NP_001093582.1).                                                                                                                                                        |
| <b>Immunogen Sequence</b> | NKDPLNETVWGLYQKSAMKT LAQLFSGAQTAELEGAGGGA KKGKKKKGSSFTVSALFRE NLNKLMTNLRSTHPPHFVRCI<br>IPNETKTPGAMEHELVHLQL RCNGVLEGIRICRKGFPSPRI LYADFKQRYKVLNASAIPEG QFIDSKKASEKLLASIDIDH<br>TQYKFGHTKVFFKAGLLGLL EEMRDDKLAQLITRTQARCR GLFARVEYQRMVERREAIFC IQYNIRSFMNKVHWPWMKL |



Immunofluorescence analysis of mouse skeletal muscle cells using MYH2 rabbit polyclonal antibody (STJ117487) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of rat skeletal muscle cells using MYH2 rabbit polyclonal antibody (STJ117487) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081