

Anti-PRG4 antibody (1264-1363) (STJ117575)

STJ117575

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10216
Gene Symbol	PRG4
Uniprot ID	PRG4_HUMAN
Immunogen	
Immunogen Region	1264-1363
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1264-1363 of human PRG4 (NP_001121180.2).
Immunogen Sequence	ERAIGPSQTHIRIQYSPAR LAYQDKGVLHNEVKVSILWR GLPNVVTSAISLPNIRKPDG YDYYAFSKDQYYNIDVPSRT ARAITTRSGQTLKVVWYNCP