

Anti-CIC antibody (1379-1608) (STJ117603)

STJ117603

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23152
Gene Symbol	CIC
Uniprot ID	CIC_HUMAN
Immunogen	
Immunogen Region	1379-1608
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1379-1608 of human CIC (NP_055940.3).
Immunogen Sequence	LATSPRAILGSYRKKRKNST DLDSAPEDPTSPKRMRRRS SCSSEPNTPKSAKCEGDIFT FDRGTGEAEDVLGELEYDKV PYSSLRRLTDQRRALVMQLF QDHGFFPSAQATAAFQARYA DIFPSKVCLQLKIREVRQKI MQAATPTEQPPGAEAPLPVP PPTGTAAAPAPTPSPAGGPD PTSPSSDSGTAQAAPPLPPP PESGPGQPGWEGAPQSPPPP PGPSTAATGR



Western blot analysis of lysates from NIH/3T3 cells, using CIC Rabbit polyclonal antibody (STJ117603) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081