

Anti-TKFC antibody (226-575) (STJ117616)
STJ117616

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26007
Gene Symbol	TKFC
Uniprot ID	TKFC_HUMAN
Immunogen	
Immunogen Region	226-575
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 226-575 of human TKFC (NP_056348.2). VRRIKMATADEIVKMLDHM TTTTNASHVPVQPGSSVMM VNNLGGLSFLELGIIADATV RSLEGRGVKARALVGTfMS
Immunogen Sequence	ALEMPGISLTLVVDEPLK LIDAETAAAWPNVAVSIT GRKRSRVAPAEPQEPDSTA AGGSASKRMALVLERVCSTL LGLEEHLNALDRAAGDGDGCG TTHSRAARAIQEWLKEGPPP ASPAQLLSKLSVLLLEKMGG SSGALYGLFLTAAQAQLKA



Western blot analysis of various lysates using TKFC Rabbit polyclonal antibody (STJ117616). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081