

## Anti-USP47 antibody (1128-1287) (STJ117656)

STJ117656

### GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IHC-P/ELISA     |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse        |

### PRODUCT PROPERTIES


|                            |                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                                            |
| <b>Clone ID</b>            |                                                                                                                                                       |
| <b>Concentration</b>       | Lot specific                                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                          |
| <b>Purification</b>        | Affinity purification                                                                                                                                 |
| <b>Dilution</b>            | WB:1:200-1:2000                                                                                                                                       |
| <b>Range</b>               | IHC-P:1:50-1:200<br>ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                                                                   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                        |

### TARGET INFORMATION

|                           |                                                                                                                                                                       |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 55031                                                                                                                                                                 |
| <b>Gene Symbol</b>        | USP47                                                                                                                                                                 |
| <b>Uniprot ID</b>         | UBP47_HUMAN                                                                                                                                                           |
| <b>Immunogen</b>          |                                                                                                                                                                       |
| <b>Immunogen Region</b>   | 1128-1287                                                                                                                                                             |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 1128-1287 of human USP47 (NP_060414.3).                                                 |
| <b>Immunogen Sequence</b> | VEKMKMSQLAVLSRRWKPS EMKLDPPFQEVVLESSSVDEL REKLSEISGIPLDDIEFAKG RGTFFCDISVLDIHQDLWDN PKVSTLNVWPLYICDDGAVI FYRDKTEELMELTDEQRNEL MKKESRLQKTGHRVTYSR KEKALKIYLDGAPNKDLTQD |



Immunohistochemistry analysis of USP47 in paraffin-embedded Human colon carcinoma using USP47 Rabbit polyclonal antibody (STJ117656) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using USP47 Rabbit polyclonal antibody (STJ117656) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081