

Anti-TRAPPC9 antibody (947-1246) (STJ117722)

STJ117722

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	83696
Gene Symbol	TRAPPC9
Uniprot ID	TPPC9_HUMAN
Immunogen	
Immunogen Region	947-1246
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 947-1246 of human TRAPPC9 (NP_113654.4).
Immunogen Sequence	SHVKTLEAVLNFKYSGGPGH TEGYYRNLSLGLHVEVEPSV FFTRVSTLPATSTRQCHLLL DVFNSTEHELTVSTRSSEAL ILHAGECQRMAIQVDKFNFE SFPESPGEKGQFANPKQLEE ERREARGLEIHSKLGICWRI PSLKRSGEASVEGLLNQLVL EHLQLAPLQWDVLVDGQPCD REAVAACQVGDPVRLEVRLT NRSPRSVGPFLTVPFQDH QNGVHNYDLHDTVSVFGSS



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081