

Anti-ARF3 antibody (1-181) (STJ117846)

STJ117846

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	377
Gene Symbol	ARF3
Uniprot ID	ARF3_HUMAN
Immunogen	
Immunogen Region	1-181
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human ARF3 (NP_001650.1).
Immunogen Sequence	MGNIFGNLLKSLIGKKEMRI LMVGLDAAGKTTILYKLLG EIVTTIPTIGFNVETVEYKN ISFTVWDVGGQDKIRPLWRH YFQNTQGLIFVDSNDRERV NEAREELMRMLAEDELRLDAV LLVFANKQDLPNAMNAAEIT DKLGLHSLRHRNWWYIQTCA TSGDGLYEGLDWLANQLKNK K