

Anti-ODF2 antibody (630-829) (STJ118078)

STJ118078

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4957
Gene Symbol	ODF2
Uniprot ID	ODFP2_HUMAN
Immunogen	
Immunogen Region	630-829
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 630-829 of human ODF2 (NP_702911.1).
Immunogen Sequence	IISDLRSRIEHQGDKLEMAR EKHQASQKENKQLSLKVDEL ERKLEATSAQNIEFLQVIAK REEAIHQSQLRLEEKTRCG TLARQLSAIEDARRQVEQT KEHALSKERAAQNKILDLET QLSRTKTELSQLRRSRDDAD RRYQSRLQDLKDRLEQSEST NRSMQNYVQLKSSYANVFG DGPYSTFLTSSPIRSRSPPA



Immunofluorescence analysis of HeLa cells using ODF2 rabbit polyclonal antibody (STJ118078) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of HeLa cells using ODF2 rabbit polyclonal antibody (STJ118078) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081