

Anti-CSF3 antibody (30-207) (STJ118106)

STJ118106

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1440
Gene Symbol	CSF3
Uniprot ID	CSF3_HUMAN
Immunogen	
Immunogen Region	30-207
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of human CSF3 (NP_000750.1).
Immunogen Sequence	ATPLGPASSLPQSFLKCLE QVRKIQGDGAALQEKLVSSEC ATYKLCHPEELVLLGHSLGI PWAPLSSCP SQALQLAGCLS QLHSGFLFYQGLLQALEGIS PELGPTLDTLQLDVADFATT IWQQMEELGMAPALQPTQGA MPAFASAFQRRAGGVLVASH LQSFLEVSYRVLRLHAQP